

## Trypsin In-gel Digestion Protocol

### Destaining

- Excise protein spot/ band and cut into ~1mm cubes. Transfer to a clean microcentrifuge tube.
- 250 ul water, 15 min. (can store in water for longer)
- Remove water.
- A): for silver staining: 50 ul destaining buffer, 20-30 min.
  - Destaining buffer: 3 mL 100 mM  $\text{Na}_2\text{S}_2\text{O}_3$   
3 mL 30 mM  $\text{K}_3\text{Fe}(\text{CN})_6$
- B): for coomassie blue and Sypro Ruby staining: Destain 3 times by adding 100 ul of 25 mM  $\text{NH}_4\text{HCO}_3$  in 50% acetonitrile, incubate for 15-30 minutes with occasional vortexing, and discard the liquid with each wash.
- Remove destaining solution.
- 500 ul water, 30 min. (store overnight in water if blue color is not disappearing)
- Add enough 100% acetonitrile to fully cover gel piece, incubate for 20-30 minutes, and discard acetonitrile.
  - Gel pieces should be white/ opaque and shrunken
- Dry gel pieces using Speed Vac.

### Reduction of Disulfide Bonds

- Add 100-200 ul of 10 mM dithiothreitol (DTT) in 100 mM  $\text{NH}_4\text{HCO}_3$  to the dried gel piece and incubate for 1 hour at 56 °C.

### Alkylation with Iodoacetamide

- Let the gel pieces cool to room temp and discard DTT soln.
- Add 200 ul of 55 mM (10 mg/ml) iodoacetamide in 100 mM  $\text{NH}_4\text{HCO}_3$  until gel pieces are covered.
- Wrap tube in aluminum foil to protect sample from light and incubate 45 minutes in dark at room temperature.
- Remove iodoacetamide/  $\text{NH}_4\text{HCO}_3$  solution and discard.
- Add 50-100 ul of 100 mM  $\text{NH}_4\text{HCO}_3$  and incubate for 10 minutes at room temperature. Remove and discard solution
- Dehydrate with about 100 ul of 100% acetonitrile. Remove and discard solution when gel pieces are white/ opaque.
- Re-swell in 100 ul of 100 mM  $\text{NH}_4\text{HCO}_3$  for 10 minutes. Remove and discard solution
- Dehydrate with about 100 ul of 100% acetonitrile. Remove and discard solution when gel pieces are white/ opaque.
- Use the Speed Vac to evaporate remaining solvents.

- Add enough Trypsin solution (Trypsin: protein ratio of 1:20-1:100, w/w) to re-swell gel pieces completely at 4 °C for 30 minutes.
  - If after 30 minutes, gel pieces are uncovered, add more buffer to cover gel pieces.
  - Trypsin should be 6 ng/ul for silver stain and 12 ng/ul for coomassie blue.
  - **Digestion buffer:** 2.5 mL 100% acetonitrile, 45.5 mL H<sub>2</sub>O, 2 mL 1M NH<sub>4</sub>HCO<sub>3</sub>,
- Digest overnight at 37 °C.
  - Digest no more than 18 hours

### **Peptide Extraction**

- Transfer supernatant into new labeled vial
- Solution A 30-50 ul, 1 hour at room temperature
- Transfer the solution to the same vial as 1
- Solution B 30-50 ul, 1 hour at room temperature
- Transfer the solution to the same vial as 1
- Add enough acetonitrile to dehydrate
- Transfer the solution to the same vial as 1
- Dry the solution using Speed Vac. (can store in 4°C once dry)
- **Solution A:** 5% Formic Acid in H<sub>2</sub>O
- **Solution B:** 5% Formic Acid in 50:50 H<sub>2</sub>O: Acetonitrile

### **At this point, ship samples to VGN**

### **Re-suspension**

- Add 8-10 ul loading solution (2.5% MeCN, 2.5% Formic Acid) and run on MS